

13000 rpm at 4°C for 10 min. The supernatant was carefully collected and stored at –20°C until further analyses. Protein concentrations were determined using a bicinchoninic acid (BCA) protein assay kit (Pierce™ BCA Protein Assay Kit, Thermo Fisher) according to the manufacturer's protocol.

### *Brain collection*

The brain of frozen bees was removed in total by opening the cranium using a scalpel and forceps, as described before (Christen et al. 2016). The tissue was collected into ice-cold 2x concentrated lysis buffer (100 mM Tris pH 7.4, 300 mM NaCl, 30 mM EDTA, 0.2% Triton; Duong et al. 2004)) with addition of protease inhibitor cocktail (Roche cOmplete 04693124001). The samples were stored on ice during collection. For tissue disruption the collected brains were ground using a small glass bar as pestle before lysing on ice for 30 min with two intermediate vortexing steps. Lysis was finally enhanced by sonication on ice by 1-3 cycles with 10 second bursts with 30 second cooling intervals (amplitude 100%, cycle 1 in a Branson Sonifier 450 W). Cellular debris was removed by centrifugation at 13000 rpm at 4°C for 10 min. The supernatant was carefully collected and stored at –20°C until further analyses. Protein concentrations were determined using a bicinchoninic acid (BCA) protein assay kit (Pierce™ BCA Protein Assay Kit, Thermo Fisher) according to the manufacturer's protocol.

### *Analysis of antibody binding by dot blot*

Dot blot analysis was run according to the protocols on the Licor home page ([https://www.licor.com/bio/applications/quantitative\\_western\\_blot/protocol.htm](https://www.licor.com/bio/applications/quantitative_western_blot/protocol.htm)) with slight modifications. BSA (1 µL of 2 mg/mL BSA in PBS) as negative control, secondary antibody solution (IRDye® 800CW Goat anti-Rabbit IgG, 1:500 in Licor blocking buffer diluted 1:1 with PBS-T) as positive control, hemolymph, fat body and brain lysate were spotted separately on a nitrocellulose membrane (Amersham™ Protran® Western blotting membranes, This article is protected by copyright. All rights reserved

0.45  $\mu$ m). After drying, the membrane was incubated in Licor Blocking buffer (927-40100) for 1h at room temperature with agitation. After addition of 0.2% Tween-20, the membrane was incubated with the polyclonal antibody against vitellogenin in a 1:7000 dilution (Davids Biotechnologie Regensburg) for 1 h at room temperature with agitation. After three washes with PBS-T (0.2% Tween 20) the membrane was incubated for 45 min at room temperature in secondary antibody solution (IRDye® 800CW Goat anti-Rabbit IgG (H + L), 1:40.000 in Licor blocking buffer diluted 1:1 with PBS-T). After 3 washes with PBS-T (0.2% Tween 20) the vitellogenin specific signals were acquired at 800 nm using the Licor Odyssey® CLx Imaging System.

#### *SDS Polyacrylamide gel electrophoresis*

SDS polyacrylamide gel electrophoresis was run according to the protocols on the Licor home page ([https://www.licor.com/bio/applications/quantitative\\_western\\_blot/protocol.htm](https://www.licor.com/bio/applications/quantitative_western_blot/protocol.htm)) with slight modifications. In brief, one-dimensional gel electrophoresis was carried out in vertical polyacrylamide gels (10.1 x 7.3 x 0.1 cm) containing 0.1% SDS with a 4% stacking gel on top of the separating gel. Samples were diluted with 4x concentrated Orange G sample buffer for Licor blots (250 mM Tris-HCl, pH 6.8, 12% SDS, 50% Glycerol, 6% 2-mercaptoethanol, 0.2% Orange G), heated for 5 min at 95°C and subjected to electrophoresis at constant voltage (140 V) for 20 min and 190 V until the dye front has run out. Two types of one-dimensional gels were run (10% and 7.5% polyacrylamide/ 0.1% SDS gels) for the separation of proteins in the range of 30–200 kDa, with the later found to be superior for separation of the two vitellogenin forms.

#### *Protein extraction and digestion*

The extracted proteins were resolved on one- dimensional sodium dodecyl sulfate polyacrylamide gel electrophoresis (SDS-PAGE) using a pre-cast gradient gel (4 - 20% polyacrylamide, This article is protected by copyright. All rights reserved